

Applicant: Ault-Riche *et al.*
Serial No. 09/910,120 Filed: July 18, 2001
Our Docket No.: 25885-1751

Sorting by pools: Decreasing pool diversities




FIG. 2

Sorting by pools: Screening large diversity libraries




FIG. 3

Searching a mutation library







FIG. 4

Making a recombinant antibody library



Spleen cells or PBLs



Expression



Antibodies

FIG. 5

Creating the master antibody library: Primer incorporation

1. mRNA purification from spleen or PBLs

2. cDNA synthesis

3. amplification

4. assembly





FIG. 6

Creating the master antibody library: Linker addition

1. mRNA purification from spleen or PBLs

2. cDNA synthesis

3. amplification

4. assembly





FIG. 7

Searching a recombinant antibody library




FIG. 8

Physical elements to include in the kits and combinations

- *Anti-tag Arrays™*

- Primer sets





- Readers


- Software

FIG. 9


Making the V_LFOR primers: Solid phase synthesis




1. Synthesize oligo on solid support



2. Add aminolink prior to cleavage



3. Couple to tosyl activated magnetic beads



4. Extended by hybridizing with DNA patch and ligating









FIG. 10

Making the V_{LFOR} primers: Overlapping hybridization



- Synthesize 4,028 different oligos:
 (26 for J_{kappa} for : 2,000 for Epitope, 2,000 for D; 2 for Common


- Assemble oligos for + and - strands of the different regions



- Ligase the assembled oligos




- 1st strand synthesis with biotinylated primer



- 2nd strand synthesis with non-biotinylated primer



- Bind to ovidin coated magnetic beads and then denature



- Purify non-biotinylated ssDNA




FIG. 11

Building the collection of antibody/tag pairs: Hybridoma screening




FIG. 12

Table 3 Primers for PCR Amplification of Human Antibody Variable Regions (V genes)

1. V gene primary PCR

A. Human VH back primers (sense)

HuVH1aBACK	5'-CAG GTG CAG CTG GTG CAG TCT GG-3'
HuVH2aBACK	5'-CAG GTC AAC TTA AGG GAG TCT GG-3'
HuVH3aBACK	5'-GAG GTG CAG CTG GTG GAG TCT GG-3'
HuVH4aBACK	5'-CAG GTG CAG CTG CAG GAG TCG GG-3'
HuVH5aBACK	5'-GAG GTG CAG CTG TTG CAG TCT GC-3'
HuVH6aBACK	5'-CAG GTA CAG CTG CAG CAG TCA GG-3'

B. Human JH forward primers (anti-sense)

HuJH1-2FOR	5'-TGA GGA GAC GGT GAC CAG GGT CCC-3'
HuJH3FOR	5'-TGA AGA GAC GGT GAC CAT TGT CCC-3'
HuJH4-5FOR	5'-TGA GGA GAC GGT GAC CAG GGT TCC-3'
HuJH6FOR	5'-TGA GGA GAC GGT GAC CGT GGT CCC-3'

C. Human V kappa back primers (sense)

HuVk1aBACK	5'-GAC ATC CAG ATG ACC CAG TCT CC-3'
HuVk2aBACK	5'-GAT GTT GTG ATG ACT CAG TCT CC-3'
HuVk3aBACK	5'-GAA ATT GTG TCG CAG TCT CC-3'
HuVk4aBACK	5'-GAC ATC GTG ATG ACC CAG TCT CC-3'
HuVk5aBACK	5'-GAA ACG ACA CTC ACG CAG TCT CC-3'
HuVk6aBACK	5'-GAA ATT GTG CTG ACT CAG TCT CC-3'

C. Human V lambda back primers (sense)

HuVλ1BACK	5'-CAG TCT GTG TTG ACG CAG CCG CC-3'
HuVλ2BACK	5'-CAG TCT GCG CTG ACT CAG CCT GC-3'
HuVλ3aBACK	5'-TCC TAT GTG CTG ACT CAG CCA CC-3'
HuVλ3bBACK	5'-TCT TCT GAG CTG ACT CAG GAC CC-3'
HuVλ4BACK	5'-CAC GTT ATA CTG ACT CAA CCG CC-3'
HuVλ5BACK	5'-CAG GCT GTG CTC ACT CAG CCG TC-3'
HuVλ6BACK	5'-AAT TTT ATG CTG ACT CAG CCC CA-3'

D. Human J kappa forward primers (anti-sense)

HuJk1FOR	5'-ACG TTT GAT TTC CAC CTT GGT CCC-3'
HuJk2FOR	5'-ACG TTT GAT CTC CAG CTT GGT CCC-3'
HuJk3FOR	5'-ACG TTT GAT ATC CAC TTT GGT CCC-3'
HuJk4FOR	5'-ACG TTT GAT CTC CAG CTT GGT CCC-3'
HuJk5FOR	5'-ACG TTT AAT CTC CAG TGT GGT CCC-3'

D. Human J. lambda forward primers (anti-sense)

HuJλ1FOR	5'-ACC TAG GAC GGT GAC CTT GGT CCC-3'
HuJλ2-3FOR	5'-ACC TAG GAC GGT CAG CTT GGT CCC-3'
HuJλ4-5FOR	5'-ACC TAA AAC GGT GAG CTG GGT CCC-3'

FIG. 13A

2. Linker fragment PCR

F. Reverse JH for scFv linker (sense)

FR4 heavy linker
 RHuJH1-2 5'-GC ACC CTG GTC ACC GTC TCC TCA GGT GG-3'
 RHuJH3 5'-GG ACA ATG GTC ACC GTC TCT TCA GGT GG-3'
 RHuJH4-5 5'-GA ACC CTG GTC ACC GTC TCC TCA GGT GG-3'
 RHuJH6 5'-GG ACC AGG GTC ACC GTC TCC TCA GGT GG-3'

F. Reverse Vκ for scFv linker (anti-sense)

FR1 light linker
 RHuVκ1aBACKFv 5'-GG AGA CTG GGT CAT CTG GAT GTC CGA TCC GCC-3'
 RHuVκ2aBACKFv 5'-GG AGA CTG AGT CAT CAC AAC ATC CGA TCC GCC-3'
 RHuVκ3aBACKFv 5'-GG AGA CTG GGT CAA CAC AAT TTC CGA TCC GCC-3'
 RHuVκ4aBACKFv 5'-GG AGA CTG GGT CAT CAC GAT TTC CGA TCC GCC-3'
 RHuVκ5aBACKFv 5'-GG AGA CTG GGT GAG TGT GGT TTC CGA TCC GCC-3'
 RHuVκ6aBACKFv 5'-GG AGA CTG AGT CAG CAC AAT TTC CGA TCC GCC-3'

F. Reverse Vλ for scFv linker (anti-sense)

FR1 light linker
 RHuVλBACK1Fv 5'-GG CGG CTG GGT CAA CAC AGA CTG CGA TCC GCC ACC GCC AGA G-3'
 RHuVλBACK2Fv 5'-GC AGG CTG AGT CAG AGC AGA CTG CGA TCC GCC ACC GCC AGA G-3'
 RHuVλBACK3aFv 5'-GG TGG CTG AGT CAG CAC ATA GGA CGA TCC GCC ACC GCC AGA G-3'
 RHuVλBACK3bFv 5'-GG GTC CTG AGT CAG CTC AGA AGA CGA TCC GCC ACC GCC AGA G-3'
 RHuVλBACK4Fv 5'-GG CGG TTG AGT CAG TAT AAC CTG CGA TCC GCC ACC GCC AGA G-3'
 RHuVλBACK5Fv 5'-GA CGG CTG AGT CAG CAC ATA GGA CGA TCC GCC ACC GCC AGA G-3'
 RHuVλBACK6Fv 5'-TG GGG CTG AGT CAG CAT AAA ATT CGA TCC GCC ACC GCC AGA G-3'

3. Pull-through primers for introduction of restriction sites*

G. Human VH back (Sfi) primers (sense)

HuVH1aBACKSfi FR1 heavy
 5'-GTC CTC CGA ACT GCG GCG CAG CCG GCG ATG GCC CAG CTG GTG CAG TCT GG-3'
 HuVH2aBACKSfi
 5'-GTC CTC CGA ACT GCG GCG CAG CCG GCG ATG GCC CAG GTC AAC TTA AGG GAG TCT GG-3'
 HuVH3aBACKSfi
 5'-GTC CTC CGA ACT GCG GCG CAG CCG GCG ATG GCC GAG GTG CAG CTG GTG GAG TCT GG-3'
 HuVH4aBACKSfi
 5'-GTC CTC CGA ACT GCG GCG CAG CCG GCG ATG GCC CAG GTG CAG CTG CAG GAG TCG GG-3'
 HuVH5aBACKSfi
 5'-GTC CTC CGA ACT GCG GCG CAG CCG GCG ATG GCC CAG GTG CAG CTG TTG CAG TCT GG-3'
 HuVH6aBACKSfi
 5'-GTC CTC CGA ACT GCG GCG CAG CCG GCG ATG GCC CAG GTA CAG CTG CAG CAG TCA GG-3'

H. Human J kappa forward (Not) primers (anti-sense)

HuJk1FORNot FR4 light
 5'-GAG TCA TTC TCG ACT TCG GCG CCG ACC TTT GAT TTC CAC CTT GGT CCC-3'
 HuJk2FORNot
 5'-GAG TCA TTC TCG ACT TCG GCG CCG ACC TTT GAT CTC CAG CTT GGT CCC-3'

H. Human J kappa forward (Not) primers (anti-sense)(continued)

HuJk3FORNot FR4 light
 5'-GAG TCA TTC TCG ACT TCG GCG CCG ACC TTT GAT ATC CAG CTT GGT CCC-3'
 HuJk4FORNot
 5'-GAG TCA TTC TCG ACT TCG GCG CCG ACC TTT GAT CTC CAC CTT GGT CCC-3'
 HuJk5FORNot
 5'-GAG TCA TTC TCG ACT TCG GCG CCG ACC TTT AAT CTC CAG TCG TGT CCC-3'

H. Human J lambda forward (Not) primers (anti-sense)



HuJl1FORNOT FR4 light
 5'-GAG TCA TTC TCG ACT TCG GCG CCG ACC TAG GAC GGT CAC CTT GGT CCC-3'
 HuJl2-3FORNOT
 5'-GAG TCA TTC TCG ACT TCG GCG CCG ACC TAG GAC GGT CAG CTT GGT CCC-3'
 HuJl4-5FORNOT
 5'-GAG TCA TTC TCG ACT TCG GCG CCG ACC TAA AAC GGT GAG CTG GGT CCC-3'

*Recognition site for restriction enzyme is underlined.



step I

Tag and assemble immunoglobulin genes



Create 1,000 sub-libraries by separate PCR amplification
 reactions using tag-specific PCR primers



FIG. 14A



step II




FIG. 14B

step III





FIG. 14C



summary



Round	Library diversity	Sub-libraries	Spot diversity	Arrays used
0	1,000,000,000	1,000		
I	1,000,000		1,000	1,000
II	1,000		1	1

FIG. 14D

Modification searches




FIG. 15

Simultaneous searches

Round Arrays Boit Probe

I 1,000 Abs Aqs

II 1,000 Abs Aqs

III $\frac{1,000 \text{ Aqs Abs}}{3,000}$






FIG. 16

Protein interaction mapping



Spots can contain mixtures of enzymes
for detection or pathway engineering

FIG. 17

Protein interaction mapping




FIG. 18




FIG. 19